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1. Allow the device to warm-up for 15 minutes
2. Perform an offset calibration of the pressure and flow sensors (**Zero!** calibration) on the warmed-up device.

- For a short time after turning on the instrument some displays may vary slightly from zero until the optimal operating temperature is reached (10 to 15 Min).
- Therefore the offset calibration (**Zero!** calibration) must not be performed, as long the instrument is cold.
- During a zero calibration it is important that no pressure or flow is applied to any connector and the yellow protection caps are not applied.

3. Check the gas settings in the menu (Gas Type/Standard).

Gas standard is set by default to ATP (Ambient temp. & pressure).
Selecting the wrong gas type or wrong gas standard can lead to measurement inaccuracies of up to 20 %

4. If you want to measure respiratory parameters, check the trigger settings in the menu (Trigger).

Trigger settings are set by default to "Adult"

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FlowLab Software FlowLab Software

View your flow, pressure and volume measurements comfortably on your computer monitor, while simultaneously saving respiratory parameters.

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FlowLab software download:
www.imtanalytics.com/FlowLab

EasyCal annual Recalibration EasyCal jährliche Rekalibrierung

EasyCal

To ensure the accuracy of your measurements, you must recalibrate your FlowAnalyser annually at IMT Analytics.

Our EasyCal recalibration service makes recalibration fast and easy! Worldwide express 10-day service including shipping. For all details visit: www.imtanalytics.com/easycal

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User Manual Gebrauchsanweisung

Please read the user manual carefully before using the FlowAnalyser.
Online at: www.imtanalytics.com/analyser_manuals

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USB-Cable
USB Kabel

Calibration certificate
Kalibrationszertifikat


Technical Support Technischer Support


IMT Analytics

IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland
T +41 81 750 67 10 . www.imtanalytics.com . technicalsupport@imtanalytics.com

Operating Elements Bedienelemente

- 1 **Low flow measurement channel**
Messkanal für niedrige Flüsse
- 2 **High flow measurement channel**
Messkanal für hohe Flüsse
- 3 **Display**
Bildschirm
- 4 **Context-sensitive access keys**
Kontextsensitive Bedientasten
- 5 **Navigation arrow keys**
Pfeiltasten zur Navigation
- 6 **Support handle**
Tragegriff
- 7 **USB port**
USB Schnittstelle
- 8 **Device fan**
Gerätelüfter
- 9 **Direct Access Control (DAC) keys**
Direct Access Control (DAC) Tasten
- 10 **PF-300:**
Differential pressure measurement port 1 (against port 2)
Differenzdruck Messanschluss 1 (gegen Anschluss 2)
- 11 **PF-301 VAC:**
+/- 1 bar pressure and vacuum measurement port
+/- 1 bar Druck und Vakuum Messanschluss
- 12 **PF-302 LOW:**
Low pressure measurement port
Niederdruck Messanschluss


- 11 **PF-300:**
Differential pressure measurement port 2 (against port 1)
Differenzdruck Messanschluss 2 (gegen Anschluss 1)
- 12 **PF-301 VAC:**
Differential pressure measurement port (against ambient pressure)
Differenzdruck Messanschluss (gegen Umgebungsdruck)
- 13 **PF-302 LOW:**
Differential pressure measurement port (against ambient pressure)
Differenzdruck Messanschluss (gegen Umgebungsdruck)
- 14 **High pressure measurement port**
Hochdruck Messanschluss
- 15 **Power indicator LED**
Versorgungs-Anzeige LED
- 16 **On/Off switch**
Ein/Aus Schalter
- 17 **RS 232 port**
RS 232 Schnittstelle
- 18 **Power supply connector**
Netzteil Anschluss
- 19 **External trigger interface**
Externer Trigger Schnittstelle

Measurement Setup Messaufbau

! In order to protect the instrument from contaminations with particles in the air it is important to use the protection filter, which is delivered with each instrument. The protection filter has to be used on the High flow as well as on the Low flow channel in upstream position in front of the device. Before use: make sure that all system components (tube, protection filter, test lungs, etc.) are clean and free of dust.

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

Use the Context-sensitive access keys for general operation.

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
Use the Direct Access Control (DAC) keys located next to each mechanical connection to display information and current measurement.

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
Use navigation arrow keys to move the cursor marker in setting screens.

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The display light contrast can be adjusted to your preference and viewing angle by pressing and holding the first button from the left then use the arrow key to select the desired contrast level.

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There are no additional sensors in the low flow channel to measure temperature, humidity, or O₂ concentration which have an impact on flow measurement. To achieve very accurate measurements it is recommended to connect the back end of the low flow channel to the high flow channel.

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