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
EasyLung Neonatal 500.040.000

Manufacturer and Technical Support

Hersteller und technischer Support
Fabricant et support technique
Fabricante y asistencia técnica
Produttore e assistenza tecnica

English – EN

Any serious incidents that have occurred in relation to the EasyLung Neonatal must be reported to the manufacturer.

Intended use

The EasyLung Neonatal is intended to be used exclusively for regular functional testing of ventilators or anaesthetic machines for neonates in manufacturing facilities, hospitals or other healthcare facilities.

If you have received the test lung with damaged packaging, please contact the technical support listed in this instruction for use.

To maintain the full functionality of the EasyLung Neonatal, a tube system with a taper of 15mm diameter should be used.

The EasyLung Neonatal shall hang on its suspension during operation.

Cleaning Process

Before reusing, the EasyLung Neonatal must be checked for visible external damage and must be tested for function. Make sure the EasyLung Neonatal is protected against mechanical damage during transport.

 Don't use a damaged EasyLung Neonatal.

Visible contamination on the surface of the test lung should be removed with a clear and lint-free cloth (use e.g. Gigasept FF from Schülke & Mayr GmbH).

The complete EasyLung Neonatal can be autoclaved according to the guidelines and recommendations for hygiene preparation at 134 °C and 2bar for a duration of 3.5min as well as a drying time of at least 2min. The EasyLung Neonatal must be cooled down at room temperature for at least 120 min.

 Don't use the EasyLung Neonatal for:
Gas sterilization, Hot air sterilization, Radiation sterilization

Contraindications

The test lung is intended exclusively for use with ventilators and anaesthetic machines for neonates. Direct use on patients is not permitted. Therefore, no contraindications are defined.

Specifications

Volume of the silicone bag:	max. 50 mL	Temperature range for operation: for transport and storage:	5 °C–40 °C – 10 °C–60 °C
Compliance:	1 mL/mbar ± 5 %	Humidity range for operation: for transport and storage:	10 % RH–95 % RH 5 % RH–95 % RH
Resistance:	2440 mbar/L ² /s ² ± 10 %	Atmospheric pressure range for operation: for transport and storage:	600 hPa–1100 hPa 540 hPa–1100 hPa
Useful life time:	1'000'000 breathing cycles or max. 50 reprocessing cycles, whatever comes first		

Free from latex, phthalates and BPA

 The EasyLung Neonatal is made of silicone. During preparation, ensure that no particles enter into the test lung and the test lung does not get damaged. IMT Analytics AG is not liable for damages caused by disregarding this instruction for use.


CE conformity


This test lung satisfies the EU Medical Device Regulation (2017/745/EU).

UKCA Conformity

This test lung satisfies the regulation of United Kingdom Conformity Assessment.

Explanation of symbols

	Manufacturer's address incl. date of manufacture
	Catalogue number
	Lot number
	MDR Importer: IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein
	Mark for United Kingdom Conformity Assessment
	CE marking
	Medical device
	Unique device identification
	Date Stamp
	Consult instructions for use

	WEEE marking Waste Electrical and Electronic Equipment
	Keep dry
	Keep away from sunlight
	European Authorized Representative: IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein
	UK Responsible Person: Mr. Friberg 104B Elers road, London, W13 9QE England United Kingdom
	Temperature limitation for transport and storage
	Humidity limitation for transport and storage
	Atmospheric pressure limitation for transport and storage
	Sterilizable in a steam sterilizer (autoclave) at 134 °C

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Importateur MDR : IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein
Représentant européen autorisé : IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein