SmartLung 2000
SmartLung Adult
SmartLung Infant

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User manual

Manual de instrucciones

Release 03, 2022-11

Deutsch	3
English	17
Español	31

Einleitung

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Preface

The SmartLung 2000, SmartLung Adult and the SmartLung Infant are the same size as conventional test lungs but offer both an adjustable compliance and a resistance that can be varied between Rp5 and Rp200. The adjustable leakage also enables the user to verify the leak compensation of the ventilator. The patient trigger functions can also be verified.

The SmartLung is the ideal tool for verifying ventilators and anaesthesia equipment and is the ideal supplement for CITREX and FlowAnalyser flow and pressure measurement devices.

Typical settings

	SmartLung 2000 2L	SmartLung Adult 1L	SmartLung Infant 0.5L
Adult			
Compliance (mL/mbar)	25, 60, 75	10, 15, 20, 30	
Compliance (mL/kPa)	200, 500		
Resistance	Rp5, Rp20, Rp50	Rp5, Rp20, Rp50	
Paediatric			
Compliance (mL/mbar)		10, 15, 20	
Resistance		Rp20, Rp50	
Infant			
Compliance (mL/mbar)			1, 2, 3, 5
Resistance			Rp50, Rp200

Technical Data

	SmartLung 2000 2L	SmartLung Adult 1L	SmartLung Infant 0.5L
Compliance (mL/mbar)	25, 60, 75	10, 15, 20, 30	1, 2, 3, 5
Resistance	Rp5, Rp20, Rp50, Rp200		
Volume (mL)	0–1000 2 L bag	0–600 1 L bag	0–200 mL 0.5 L bag
Leak (L/min)	0–10		
Weight (g)	400	325	285
Dimensions (mm)	400×140×40	303×115×40	270×115×40
Deadspace (mL)	3.7		
Components	replaceable		
Cleaning	by disinfection		

* Tolerances: Resistance \pm 10%, compliance \pm 20%

Approvals

The test lungs meet CE guidelines.

CE conformity

These test lungs satisfies the EU Medical Device Regulation (2017/745/EU).

UKCA Conformity

These test lungs satisfies the regulation of United Kingdom Conformity Assessment.

Resistance curves


Due to the non-laminar flow, the resistance R_p (flow resistance, parabolic resistance) depends on the flow. As the the flow rate increases, the resistance increases linearly and the pressure increases as the square of the flow rate. It therefore behaves similarly to the resistance of an endotracheal tube.

$$\text{Pressure drop} \quad P = K \times \text{flow}^2$$


$$\text{Resistance} \quad R = K \times \text{flow}$$

With K defined for each R_p :


	R_{p5}	R_{p20}	R_{p50}	R_{p200}
K	2.7	17.61	108.7	2440

The pressure drop P for various R_p values at different flow rates (parabolic curve)


Resistance with different flow rate (R_p5 and R_p20)


Resistance with different flow rate (R_p50 and R_p200)

Compliance curves


The compliance (stretchability) of the SmartLung test lungs changes for different tidal volumes. The specified compliance values therefore apply to the particular volume indicated.

Compliance curve for SmartLung Infant 0.5L

Compliance curve for SmartLung Adult 1L

Compliance curve for SmartLung 2000 2L

Packaging and transport

The SmartLung is supplied in a carry bag that can be also be used to transport a maximum of two lungs.

Adjusting the compliance

Using the slider, various lung compliance values can be set without using an adapter. To make the adjustment, press the lobes in the middle of the lung together and push the slider into the desired position with the other hand.

Adjusting the resistance

Different resistances can be simulated by simply turning the connector.

Adjusting the leakage

The leakage can be steplessly adjusted by turning the setting screw.

Fixing

Next to the pins for the bag there is a hole in the lung lobes that is used to fix the test lung.

Using the SmartLung for tests in accordance with the standard

Ventilators and anaesthesia equipment must be verified in accordance with the following standards:

- EN 60601-2-12 Critical care ventilators (substitute for DIN EN 794-1:2001-02)
- EN 60601-2-13 Anaesthetic systems (replacement for DIN EN 740:1999-02)
- EN 794-3 Transport ventilators
- ASTM F1100-90 (1997) Critical care ventilators (US FDA)

These standards determine the following settings for the SmartLung:

	Adult	Paediatric	Infant, neonatal
Tidal volume Vt (mL)	500	300	30
Rate (bpm)	10	20	30
I:E ratio	1:2	1:2	1:2
Test lung used	SmartLung 2000 2L	SmartLung 2000 2L	SmartLung Infant 0.5L
Compliance (mL/mbar)	50	20	1
Resistance Rp	5	20	50

SmartLung settings for tests in accordance with the standard

Scope of delivery/Replacemen parts

	SmartLung 2000 2L	SmartLung Adult 1L	SmartLung Infant 0.5L
Complete	302.730.100	300.162.100	300.400.104
Carry bag	302.733.000	300.360.000	
Silicone bag	302.918.000	301.185.000	300.764.000
Fixing pins for the silicone bag	302.915.000	300.317.000	
Compliance slider	302.914.000	300.300.000	
Resistance head with Lobes	304.043.000	-	
Resistance head	-	300.314.000	
Lobes	-	700.072.000	700.069.000
Lock bar for the lobes	-	300.318.000	
User manual	700.068.000		

SmartLung 2000


SmartLung Adult and SmartLung Infant

Disassembly/assembly of the SmartLung 2000

Remove fixing pins for the silicone bag.

Remove lock bar.

Detach blue bag head from the resistance head.

Pull out silicone bag.

The lung is assembled in the reverse order.

Disassembly/assembly SmartLung Adult SmartLung Infant

Remove the lock bar.

Push the resistance head downwards.

Pull the silicone bag out backwards and remove the pins.

Pull the resistance head out sideways and remove the whole unit.

The silicone bag can be detached from the resistance head.

The lung is assembled in the reverse order. Finally, push the resistance head forwards and insert the lock bar.

Cleaning

All parts of the lung can be cleaned with alcohol-based disinfectants (e. g. Pantasept Spray from Xeropharm).

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IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland