

EasyLung Neonatal

Test lung simple

testlung

the unique simulation

The EasyLung Neonatal is the ideal tool to test neonatal and infant ventilators.



Technical Specification

Parameter	Spezifikation
Volume	50mL
Compliance	1 mL/mbar @Vt 20mL
Resistance	Rp 200
Dimensions	45mm × 150mm
Weight	25g
Reprocessing	Autoclavable

IMT.Analytics

IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland
T +41 81 750 67 10 . www.imtanalytics.com